

(19) World Intellectual Property
Organization
International Bureau



(43) International Publication Date
6 May 2005 (06.05.2005)

PCT

(10) International Publication Number
WO 2005/041148 A1

(51) International Patent Classification⁷: **G09B 23/28**
(21) International Application Number:
PCT/US2003/030128

(74) Agent: **BOOKS, Mark, E.**; Polster, Lieder, Woodruff & Lucchesi, L.C., 12412 Powerscourt Drive, Suite 200, St. Louis, MO 63131-3615 (US).

(22) International Filing Date:
25 September 2003 (25.09.2003)

(25) Filing Language: English

(26) Publication Language: English

(71) Applicant (for all designated States except US): **EVER-EST BIOMEDICAL INSTRUMENTS [US/US]**; 16690 Swingley Ridge Road, Chesterfield, MO 63017 (US).

(81) Designated States (*national*): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NI, NO, NZ, OM, PG, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SY, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.

(84) Designated States (*regional*): ARIPO patent (GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

(72) Inventors; and


(75) Inventors/Applicants (for US only): **CAUSEVIC, Elvir [US/US]**; 16315 Autumn View Terrace, Ellisville, MO 63011 (US). **POPELKA, Gerald, R. [US/US]**; 21 Berkley Lane, Ladue, MO 63124 (US). **KROHN, Randall, J. [US/US]**; 16534 Lancaster Estates Drive, Wildwood, MO 63040 (US).

Published:

— with international search report

[Continued on next page]

(54) Title: HUMAN BIOELECTRIC SIGNAL SIMULATOR



(57) Abstract: A medical training device (10) provided multiple sensor contact points (20) for delivering simulated OAE and bioelectric signals to a testing or monitoring device applied to the device by a medical technician. The sensor contact points (20) in the medical training device (10) are operatively coupled to one or more signal generators and a computer system (32), whereby simulated OAE and bioelectric signals representative of selectively controlled normal and abnormal patient conditions may be provided to the applied testing or monitoring device, through the sensor contact points (20), permitting the medical technician to practice utilizing the testing or monitoring device and interpreting the results, thereby providing an improved false positive rate in practical applications.



For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.